
MEDI-PUMP[®]

 **Warning:** Thomas compressors are precision-made, and carefully assembled and wired. Therefore, do not disassemble or attempt to repair these products. Only qualified personnel should perform repair service. 

**MODEL 1205, 1206 & 1207
PORTABLE ASPIRATOR**

DEVICE SPECIFICATIONS

PUMP

12V DC oil-less piston type

PERFORMANCE

Vacuum Range: 0 to 550 mm Hg
Free Air Flow: 11 LPM

CONTROLS

Vacuum Regulator: Rotary type on panel
Vacuum Gauge: 0-30"/760 mm Hg on panel

BATTERY

Type: Rechargeable sealed lead-acid
Capacity: 12 volts, 3 ampere-hour
Average Run Time: 30 minutes when fully charged
Charge Time: 12 hours or less to 85% charging, depending on initial charge of battery.
Charge Method: A) 110 or 220 outlet with supplied charger.
B) 12 V DC vehicle outlet with optional cigarette lighter cord.

COLLECTION DEVICE

Jar: Autoclavable
Jar Lid: Autoclavable with mechanical FULL LEVEL shut-off valve
Capacity: 1200cc standard

PATIENT SIDE TUBING

Tubing: 6' Disposable plastic

INSTRUCTIONS FOR USE

IMPORTANT: This device is to be used by persons properly trained in medical suction techniques and in the operation of suction equipment. Improper use could cause injury. Thoroughly read this operation manual and familiarize yourself with the operation of the Medi-Pump before using this device.

 DANGER: POSSIBLE EXPLOSION HAZARD IF USED IN THE PRESENCE OF FLAMMABLE ANESTHETICS.

NOTE: The Medi-Pump is designed to operate with a fully charged battery (SEE CHARGING INSTRUCTIONS) or 12 V DC power cord plugged into an automobile cigarette lighter socket.

The AC wall outlet charger is for charging the battery only.

 CAUTION: DO NOT RUN UNIT WHILE CHARGER IS PLUGGED IN.


PRIOR TO INITIAL USE

Upon receiving your new Medi-Pump, perform the following initial tests to ensure that your unit is in good working order and that no damage has occurred during shipment.

1. Visually inspect all components for physical damage that may have occurred during shipping.
2. Turn the Medi-Pump ON by switching the POWER button to ON (-). Listen to verify that the vacuum pump starts. If the pump fails to start or runs slowly, charge the battery - see CHARGING INSTRUCTIONS.

3. Check the vacuum on the Medi-Pump. Turn the ON/OFF switch ON (-). Place your finger over the patient port of the collection jar and adjust the vacuum level by turning the VACUUM REGULATOR knob. Verify that the vacuum gauge reflects a change in vacuum level while turning the vacuum regulator knob. Also verify that you feel vacuum at your finger tip. Adjust the vacuum level to maximum vacuum and verify that, at maximum, the gauge indicates at least 550mm Hg. (at sea level). Turn the unit OFF by switching the POWER button to OFF (o).

4. Familiarize yourself with the location of these features: Air inlet filter, AC cord outlet, charging LED light, vacuum gauge, vacuum control knob, hydrophobic filter, vacuum and patient port of jar lid, and mechanical shut-off valve under jar lid.

OPERATING INSTRUCTIONS

1. Ensure that the vacuum source tubing, coming from the pump, is connected to the vacuum port of the collection jar. Tubing should also connect the patient port of the collection jar lid to any attachment that will be used.
2. Turn the unit ON by switching the ON/OFF switch to the ON (-) position.
3. Pinch the patient tubing closed. Using the VACUUM REGULATOR knob, adjust the vacuum level to the least amount of vacuum that will be necessary to properly suction the patient. As the VACUUM REGULATOR knob is turned clockwise - vacuum will decrease, turning it counterclockwise will increase vacuum.
4. Release the pinched tubing and suction the patient as needed.
5. The mechanical shut-off valve under the jar lid will stop suction if collected fluids reach the valve. If this happens, turn the Medi-Pump OFF (o), allow the vacuum to return to zero (0), remove the tubing from the vacuum port on the lid, lift the jar out of the support bracket, remove the lid and empty the jar. It is recommended that the jar be emptied when the fluid reaches "SAFE FULL" on the jar and after every use.
6. The hydrophobic filter is designed to filter out bacteria. It is also highly sensitive to liquids and will stop air flow if contact with liquid is made. If this happens the filter must be replaced. DO NOT operate this device without a hydrophobic filter. If this device is used without a hydrophobic filter the manufacturer is not responsible for possible damage or contamination to the pump. It is recommended that replacement parts such as hydrophobic filter, tubing and jar/lid assemblies be kept on hand at all times. See the REPLACEMENT PARTS section.

CHARGING INSTRUCTIONS

The Medi-Pump is supplied with a battery charger. To charge the battery, simply plug the AC power cord into an electrical outlet and connect the other end to the AC cord outlet on the control panel of the Medi-Pump.

NOTE:

- When charger is plugged in the LED on the control panel will be illuminated.
- Initial charging may take up to 12 hours, depending upon the condition of the battery.


CAUTION: Do not run unit while charger is plugged in.

Battery Life: A number of factors may affect the life of the battery, such as the frequency of use in terms of repeated charging and discharging, and mechanical abuse such as repeated and vigorous jarring.

It is recommended to replace the battery every 1-3 years, depending on the frequency of use. Inability to retain charge is the most detectable sign of a deteriorated battery. It is time to replace the battery when it loses its charge rapidly after a full recharge cycle. Return the Medi-Pump to your authorized service center for battery replacement.

REPLACEMENT PARTS

Component parts have been specifically designed to work with this device. When ordering replacement parts be sure to order parts designed for this device, do not substitute generic parts. The following are the parts available:

NOTE: The listed items are expendable, non-durable parts which deteriorate with time and/or exposure to cleaning solutions. These parts need to be routinely inspected and replaced by the user to keep this device in good working order, or to keep sanitary conditions. It is recommended, and the users responsibility to keep a supply of spare parts on hand.

- 638274 115V, 60Hz Charger
- 638579 230V, 60Hz Charger
- 638275 220V, 50Hz Charger

- 638507 Gasket - Jar (Not Shown)
- 638577 Carrying/Storage Case (Not Shown)
- 638163 12 Volt DC Cigarette Lighter Cord (Not Shown)

Note: This devices internal component must be serviced at an Authorized Service Center.

ROUTINE INSPECTION

To ensure that your Medi-Pump will be ready to use when it is needed, it is a good practice to check the unit out on a regular basis, or as your own Standard Operating Procedures may require. The inspection should include a few simple steps:

1. Charge the battery if necessary.
2. Check the vacuum on the Medi-Pump by placing your finger over the suction tubing or by pinching it closed, and adjust the vacuum level by turning the VACUUM REGULATOR knob. Verify that the vacuum gauge reflects a change in vacuum level while turning the vacuum regulator knob. Also verify that you feel vacuum at your finger tip. Adjust the vacuum level to maximum vacuum. Verify that, at maximum, the gauge indicates at least 550mm Hg. (at sea level). Turn the unit off by switching the "ON/OFF" switch to OFF (o).
3. Make sure that the collection jar, as well as the suction tubing and all other parts, are in good condition and properly connected.

CLEANING AND STERILIZATION FOR JAR, LID AND OTHER COLLECTION PARTS

It is recommended that the jar, lid assembly and tubing be cleaned after every use. To ensure the longest life for parts the following procedure must be followed:

1. Clean the jar, lid assembly, gasket, seal gasket, float guide and float thoroughly with mild soap. Rinse thoroughly with clean water and dry thoroughly.
2. Place parts in autoclave as follows, ensuring parts are not touching.
 - Jar with open end down.
 - Float and float guide with large ends down.
 - Gasket and seal gasket can be placed in a gusset pouch prior to autoclaving.
3. Autoclave for 3 to 5 minute cycle at 270° degrees F (132° degrees C).

NOTE: Tubing is not autoclavable and should be rinsed with clean hot tap water.

Instead of autoclave sterilization, the above parts can be soaked in an activated dialdehyde solution.

The air inlet filter in the base of the Medi-Pump should be cleaned monthly. To remove air filters, snap the round filter cover out from the bottom of the unit. The white felt filter cannot be washed. It is recommended that it be replaced every 6 months. Clean the black filter in warm mild soapy water, rinse clean and allow to air dry.

TROUBLESHOOTING

CAUTION: NEVER DISASSEMBLE THE MEDI-PUMP WHEN EITHER OF THE POWER CORDS ARE CONNECTED TO AN EXTERNAL POWER SOURCE.

Symptom #1: PUMP DOES NOT TURN ON WHEN THE POWER SWITCH IS DEPRESSED.

POSSIBLE CAUSES: 1) Battery is not charged 2) Faulty Wiring 3) Pump is not operational.

POSSIBLE SOLUTIONS: 1) Recharge the battery. 2) Make sure that the Medi-Pump is disconnected from all external power sources. Check all the electrical connections to make sure they are not loose. If all electrical connections are OK, the pump motor has probably failed. 3) Replace the pump by returning your Medi-Pump to an authorized service center.

Symptom #2: PUMP RUNS WHEN THE SWITCH IS DEPRESSED, BUT THERE IS NO VACUUM AT OUTLET.

POSSIBLE CAUSES: 1) Lid is not tight on jar. 2) Faulty tubing connections, 3) Kinked or obstructed tubing, 4) The regulator or pump is dirty, 5) Pump is not working properly.

POSSIBLE SOLUTIONS: 1) Make sure that the collection jar lid is properly seated on the jar and that the safety overflow in the lid is not stuck in the UP position. 2) Verify that all tubing connections are secure. 3) Visually inspect all tubing for obstructions or kinks. Clean or replace tubing as required. 4) and 5) Replace the pump by returning your Medi-Pump to an authorized service center..

Symptom #3: VACUUM GAUGE READS ZERO WHEN VACUUM IS AVAILABLE AT THE VACUUM OUTLET.

POSSIBLE CAUSES: 1) Faulty tubing connections 2) Gauge is broken.

POSSIBLE SOLUTIONS: 1) If tubing connections are verified to be leak free and the vacuum gauge still reads zero, it is most likely that the gauge is broken. 2) Replace the gauge by returning you Medi-Pump to an authorized service center.

Symptom #4: BATTERY WILL NOT TAKE A CHARGE.

POSSIBLE CAUSES: 1) Battery charger not working 2) Battery is old and cannot retain charge.

POSSIBLE SOLUTIONS: 1) Replace the battery charger 2) Replace battery by returning your Medi-Pump to an authorized service center.

Symptom #5: UNIT WILL NOT OPERATE ON DC POWER THROUGH THE 12V DC POWER CORD.

POSSIBLE CAUSES: 1) 12 VDC power cord improperly connected, 2) Fuse in power cord has failed.

POSSIBLE SOLUTIONS: 1) Check connection of cord to Medi-Pump 2) Replace fuse in cigarette lighter plug.

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